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Pyramids QC Final Inspection Standard

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GCX Mounting Solutions
VHM-P (Non-Locking) and VHM-PL (Locking) Variable Height Arm with Side-In Mounting Plate
(Refer to qualified personnel)

Caution: The weight of the device being mounted does not exceed the load rating of the Arm. Check the bottom surface of the Arm for the "Load Range" label.

Counter Balance Load Ranges:
25 - 45 lbs / 11.3 - 20.4 kg - Maximum Torque Weight Limit: 38 lbs / 17.2 kg
or
15 - 30 lbs / 6.8 - 13.6 kg

Installation Time:
Approximately 10-15 minutes

The purpose of this manual is to describe general installation, operation, and adjustment procedures for VHM-P and VHM-PL Series Arms. This manual should be used in conjunction with any instrument specific installation guides. Please read this manual and all instrument specific installation material before installing or using the product.

Table of Contents:

1.8 Parts Reference	2
2.8 Tools Required	2
3.8 Product Warning - Safety Labels	3
4.8 Installation, Operation and Maintenance Warnings	4
5.8 Installing the Arm in Channel	5
6.8 Mounting the Device	6
7.8 Adjusting the Arm	7
8.8 Cable Management	8
9.8 Installing Pin Cover	9
10.8 Channel Cover Installation	10
11.8 Routine Maintenance	11
12.8 Cleaning the Arm	12
13.8 Troubleshooting the VHM-P and VHM-PL	13

DU-WS-0012-01 Rev B 5/24/16 GCX Corp. Page 1 of 15

GCX Mounting Solutions
VHM-P (non-locking) and VHM-PL (locking) Variable Height Arm Operation Guide

Safety Label Locations (VHM-P Series Arms)

Refer to Installation Guide for additional information.

Remove the Device only when the Arm is at the highest position.

Do not remove device when arm is in a lowered position.

Potential pinch point that may cause personal injury.

Maximum weight for mounting plate with all function:

Counter Balance Load Range:
25 - 45 lbs / 11.3 - 20.4 kg
or
15 - 30 lbs / 6.8 - 13.6 kg

Safety Label Locations (VHM-PL Series Arms)

Refer to Installation Guide for additional information.

Potential pinch point that may cause personal injury.

Maximum weight for mounting plate with all function:

Counter Balance Load Range:
25 - 45 lbs / 11.3 - 20.4 kg
or
15 - 30 lbs / 6.8 - 13.6 kg

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